


FOR USE ON THE JAY GS[®] CUSHION

1. Remove the cover and fluid pad from the cushion to determine correct placement of the Hip Guide.
2. Detach extra Velcro[®] strips from the Hip Guide. Discard the 2 x 6" loop piece and the 2 x 6-1/4" loop piece. These pieces are not required for use with the JAY GS Cushion.
3. Place Hip Guide on cushion base at desired location. Attach the large Velcro[®] Suspender (4 x 6") so that the suspender contacts both the 2 x 6" piece of loop Velcro on the Hip Guide surface and the loop Velcro on the cushion base. Press down firmly on the suspender.

FOR USE ON THE JAY[®] CUSHION

1. Remove the cover and fluid pad from the cushion to determine correct placement of the Hip Guide.
2. Place Hip Guide at desired position on Jay Cushion base. Locate the loop Velcro on bottom of Hip Guide where it will contact the cushion base. Remove the backing to the hook and loop Velcro and press Hip Guide back into position (see diagram).
3. With Hip Guide in place, take the extra 2 x 6" loop Velcro, peel off paper backing and attach next to the Hip Guide on the base. Then take the Velcro Suspender (4 x 6") and attach it to the loop on the Hip Guide and loop on the base. Press down firmly to hold the Hip Guide in position.



Sunrise Medical • 2842 Business Park Ave. • Fresno, CA 93727 • USA
In Canada (800) 263-3390