INSTALLATION AND ADJUSTMENT OF THE SWING-AWAY LATERAL TRUNK SUPPORT

The Swing-away Lateral Trunk Support can be mounted to the lateral back plates using 2x T-nuts, 2x washers, and 2x M6x1x12 bolts per side.


- I. Location and Attachment 2
 - a. Using a 4mm hex key from the rear of the seating shell, loosen the screws(A) for the lateral back mount plate (B).
 - b. Slide the Mount plate to a position that allows some open space in the horizontal slots, so that you can attach the swing-away trunk supports.
- 2. Attaching and adjusting the Lateral Trunk Support 2 3
 - Using the 4mm hex key, mount the Lateral Trunk Support bracket holding the T-nuts in place from the rear of the backshell, thread the M6 screws (C) thru the Lateral bracket and into the T-nuts.
 - Adjust the bracket and pad mount(D) to it's optimum support position, and tighten the hardware.
 - c. Mount the pad to the pad mount(D).
 - d. Repeat the steps on both sides.
- 3. Setting Angle for Swing-away Trunk Support 4 5 6 7 8
 - a. Unlock the Swing-away lateral and swing it out of the way to gain access to the lock-screws.
 - b. Using the 4mm hex key, loosen the button head lock screws.
 - Adjust the bracket and pad to it's optimum support position, and unlock the Lateral mechanism.
 - d. Hold the adjustment in place while tightening the lock-screws.
 - e. Repeat the steps on both sides.





UNLOCK




LOOSEN









ADJUST UNLOCK AGAIN



TIGHTEN

7