

Divided Footrest Extension Mount Ordering Guide

Z-finity Hangers Multi Angle Adjustable Footplates w/ Ext & HL (Divided Footrest)								
(Order 1 each)								
Footrest width	Short Footplate Clamp (2 Hole)				Long Footplate Clamp (3 Hole)			
	Short Hemi	Short	Medium	Long	Short Hemi	Short	Medium	Long
10"	136189	136207	136225	136243	136261	136279	136296	136314
	136190	136208	136226	136244	136262	136280	136297	136315
11"	136191	136209	136227	136245	136263	136281	136298	136316
	136192	136210	136228	136246	136264	136282	136299	136317
12"	136193	136211	136229	136247	136265	136283	136300	136318
	136194	136212	136230	136248	136266	163269	136301	136319
13"	136195	136213	136231	136249	136267	136284	136302	136320
	136196	136214	136232	136250	136268	136285	136303	136321
14"	136197	136215	136233	136251	136269	136286	136304	136322
	136198	136216	136234	136252	136270	136287	136305	136323
15"	136199	136217	136235	136253	136271	136288	136306	136324
	136200	136218	136236	136254	136272	136289	136307	136325
16"	136201	136219	136237	136255	136273	136290	136308	136326
	136202	136220	136238	136256	136274	136291	136309	136327
*17"	136203	136221	136239	136257	136275	136292	136310	136328
	136204	136222	136240	136258	136276	136293	136311	136329
*18"	136205	136223	136241	136259	136277	136294	136312	136330
	136206	136224	136242	136260	136278	136295	136313	136331

***NOTE: 17" and 18" width X'cape is ONLY achieved when the dealer reattaches the x-tube assembly to the second hole on the tabs of one side; frame tabs for the 17" wide and for the 18" wide do the same thing on the other side frame.**