

FRONT ANTI-TIP REFERENCE GUIDE

(04/2008)

CASTER	3 1/4" FORK	3 1/4"+ 3/4" FORK and STEM	3 1/4"+ 1 1/2" FORK and STEM
3" MICRO	880488 880489	880488 880489	880488 880489

CASTER	4 1/4" FORK	4 1/4"+ 3/4" FORK and STEM	4 1/4"+ 1 1/2" FORK and STEM
3" MICRO	880488 880489	880488 880489	880488 880489
5" REGULAR	880488 880489	880488 880489	880490 880491

CASTER	6" FORK	6"+ 3/4" FORK and STEM	6"+ 1 1/2" FORK and STEM
3" MICRO	880490 880491	880490 880491	880494 880495

5" REGULAR	880494 880495	880494 880495	880494 880495
6" REGULAR	880494 880495	880494 880495	880494 880495
8" REGULAR	880494 880495	880494 880495	880498 880499

CASTER	5" FORK	5"+ 3/4" FORK and STEM	5"+ 1 1/2" FORK and STEM
3" MICRO	880488 880489	880490 880491	880490 880491
5" REGULAR	880490 880491	880490 880491	880494 880495
6" REGULAR	880490 880491	880490 880491	880494 880495
8" REGULAR	880494 880495	880494 880495	880494 880495

CASTER	7" FORK	7"+ 3/4" FORK and STEM	7"+ 1 1/2" FORK and STEM
3" MICRO	880494 880495	880494 880495	880494 880495
5" REGULAR	880494 880495	880494 880495	880498 880499
6" REGULAR	880494 880495	880494 880495	880498 880499
8" REGULAR	880498 880499	880498 880499	SPECIALIZED SPECIALIZED

CASTER	8" x 2" FORK	8" x 2 + 3/4" FORK and STEM	8" x 2 + 1 1/2" FORK and STEM
8" x 2"	880494 880495	880494 880495	880498 880499

KIDZ REVERSE CONFIGURATION			
CASTER	8" x 2" FORK	8" x 2 + 3/4" FORK and STEM	8" x 2 + 1 1/2" FORK and STEM
8" x 2"	880488 880489	880488 880489	880488 880489

Remarks