



Inset formula:

0" inset: (seat width) = (footrest width)

1" inset: (seat width) - 2" = (footrest width)

2" inset: (seat width) - 4" = (footrest width)

Tubular w/ Flared Cover w/o leg strap

(Order 1 each)

Width	Short/ Standard	Long
11"	050270-001	050270-002
12"	050270-003	050270-004
13"	050270-005	050270-006
14"	050270-007	050270-008
15"	050270-009	050270-010
16"	050270-011	050270-012
17"	050270-013	050270-014
18"	050270-015	050270-016
19"	050270-017	050270-018
20"	050270-019	050270-020