














Wheel Type	Axle Type		Axle P/N	Spacer P/N
Mag Wheels	Quick Release	<div style="border: 1px solid black; padding: 2px; display: inline-block;">1A</div> 	922356 Includes Spacer	<div style="border: 1px solid black; padding: 2px; display: inline-block;">2A</div>  920117
	Quad Quick Release	<div style="border: 1px solid black; padding: 2px; display: inline-block;">1B</div> 	920171	<div style="border: 1px solid black; padding: 2px; display: inline-block;">2A</div>  920117
Standard Spoke Wheels & Spinergy Wheels	Quick Release	<div style="border: 1px solid black; padding: 2px; display: inline-block;">1C</div> 	102064	N/A
	Quad Quick Release (N/A with Spinergy)	<div style="border: 1px solid black; padding: 2px; display: inline-block;">1D</div> 	922170	N/A
Quickie Performance Wheels	Quick Release (Steel)	<div style="border: 1px solid black; padding: 2px; display: inline-block;">1E</div> 	920130	N/A
	Quad Quick Release (Steel)	<div style="border: 1px solid black; padding: 2px; display: inline-block;">1F</div> 	920174	N/A
	Quick Release (Ti)	<div style="border: 1px solid black; padding: 2px; display: inline-block;">1G</div> 	920137	N/A
	Quad Quick Release (Ti)	<div style="border: 1px solid black; padding: 2px; display: inline-block;">1H</div> 	920173	N/A
Ultra Light Wheels (N/A with Quad Release)	Quick Release	<div style="border: 1px solid black; padding: 2px; display: inline-block;">1K</div> 	102064 Order Spacer	<div style="border: 1px solid black; padding: 2px; display: inline-block;">2C</div>  103093 <div style="border: 1px solid black; padding: 2px; display: inline-block;">2B</div>  100768-001

Pos.	Item Number	Description	Remarks
1A	922356	QR AXLE BIG BUTTON W/SPACER	
1B	920171	QR AXLE MAG W/QUAD NUT	
1C	102064	QR AXLE M12 NUT CLR BIG BTTN	
1D	922170	QR AXLE W/QUAD NUT PUSH BUTTON	
1E	920130	QR AXLE NON-RECESS HUB W/NUT	
1F	920174	QR AXLE NON RECESS W/QUAD NUT	
1G	920137	QR AXLE W/NUT TI	
1H	920173	QR AXLE TI W/QUAD NUT	
1I	102064	QR AXLE M12 NUT CLR BIG BTTN	
2A	920117	MAG SPACER 6S W/SET SCREW	
2B	100768-003	SPACER AXLE 220DX12.71X12MM	
2C	103093	O-RING 1/2"IDX11/16"OD BUNA-N	