



DESCRIPTION	DIM "B"
STD	7 3/8
LONG	8 7/8

Inset formula:

0" inset: (seat width) = (footrest width)

1" inset: (seat width) - 2" = (footrest width)

2" inset: (seat width) - 4" = (footrest width)

**Platform Footplate w/o leg strap**

**NOTE:** uses new style leg strap

Width	Short/ Standard	Long
11"	050249-001 (1ea) * 500690 (1ea)	050250-001 (1ea) * 500690 (1ea)
12"	050249-001(1ea) * 500690 (1ea)	050250-001(1ea) * 500690 (1ea)
13"	050249-001 (1ea) * 500691 (1ea)	050250-001 (1ea) * 500691 (1ea)
14"	050249-002(1ea) * 500692 (1ea)	050250-002 (1ea) * 500692 (1ea)
15"	050249-002 (1ea) * 500692 (1ea)	050250-002 (1ea) * 500692 (1ea)
16"	050249-003 (1ea) * 500692 (1ea)	050250-003(1ea) * 500692 (1ea)
17"	050249-003 (1ea) * 500693 (1ea)	050250-003 (1ea) * 500693(1ea)
18"	050249-004 (1ea) * 500693 (1ea)	050250-004 (1ea) * 500693 (1ea)
19"	050249-004 (1ea) * 500693(1ea)	050250-004 (1ea) * 500693(1ea)
20"	050249-004 (1ea) * 500693(1ea)	050250-004 (1ea) * 500693 (1ea)

**NOTE:** The footplate is sold separately from the Footplate Assembly.