



Inset formula:

0" inset: (seat width) = (footrest width)

1" inset: (seat width) - 2" = (footrest width)

2" inset: (seat width) - 4" = (footrest width)

Tubular w/ Flat Cover w/o leg strap (Order 1 each)		
Width	Short/ Standard	Long
11"	050024-001	050024-002
12"	050024-003	050024-004
13"	050024-005	050024-006
14"	050024-007	050024-008
15"	050024-009	050024-010
16"	050024-011	050024-012
17"	050024-013	050024-014
18"	050024-015	050024-016
19"	050024-017	050024-018
20"	050024-019	050024-020