


## Jay Xtreme Cushion Measuring Guide



PART NUMBER	CHAIR SIZE	ACTUAL WIDTH	ACTUAL LENGTH	FEMORAL SUPPORT LENGTH	PELVIC CONTOUR LENGTH	PELVIC CONTOUR WIDTH	OVERALL HEIGHT	PELVIC CONTOUR DEPTH
		A	B	C	D	E	F	G
944	14 x 14	13.75	15	4.75	5.50	12.00	2.13	1.25
946	14 x 16	13.75	17	7.00	5.50	12.00	2.13	1.25
955	15 x 15	14.75	16	6.25	5.50	13.75	2.25	1.25
956	15 x 16	14.75	17	7.25	5.50	13.50	2.25	1.25
966	16 x 16	15.75	17	7.75	4.69	13.50	2.13	1.25
967	16 x 17	15.75	18	6.75	6.75	13.25	2.25	1.25
968	16 x 18	15.75	19	8.00	6.50	13.50	2.25	1.25
976	17 x 16	16.75	17	6.50	6.25	14.50	2.25	1.25
977	17 x 17	16.75	16	7.00	6.50	14.25	2.25	1.25
986	18 x 16	17.75	17	6.50	6.25	14.75	2.25	1.25
988	18 x 18	17.75	19	8.25	6.25	13.50	2.25	1.25
960	16 x 20	15.75	21	10.25	7.00	13.50	2.25	1.25
980	18 x 20	17.75	21	10.50	7.00	13.50	2.25	1.25
900	20 x 20	19.75	21	9.50	8.00	16.50	2.25	1.25

All listed measurements are in inches