

## Caster / Drive Wheel Combinations for Pulse / ZM310

\*Original Pulse S/N PLS that did not describe frame color in the order detail were Silver Frame

\*\*Caster hub color correlates to frame color.

So if Silver Frame -hubs were grey / Black Frame - hubs were black

\*\*\*7" casters Front / 6" Casters Rear

### CASTER & FORK INFORMATION

Caster / Fork Part #	Description	Fork Color	Tire Color	Hub Color	Size	CHAIR LOCATION
116173	CASTER FORK ASSM FRNT 7IN	Black	Grey	Black	7	FRONT
127207	CASTER FORK ASSM FRNT	Black	Black	Black	7	FRONT
127206	CASTER FORK ASSM RR	Black	Black	Black	6	REAR
116175	CASTER FORK ASSM RR	Black	Grey	Black	6	REAR
126536	CASTER WHEEL 7 IN BLACK	n/a	Black	Black	7	FRONT
101691	CASTER WHEEL	n/a	Grey	Black	7	FRONT
126535	REAR CASTER WHEEL	n/a	Black	Black	6	REAR
102797	CASTER 6" RHAPSODY	n/a	Grey	Black	6	REAR
104458	REAR CASTER WHEEL GREY	n/a	Grey	Grey	6	REAR
104457	FRONT CASTER WHEEL GRAY	n/a	Grey	Grey	7	FRONT
104443	CSTR FORK ASSM FRNT S 7IN SLVR	Silver	Grey	Grey	7	FRONT
104444	CSTR FORK ASSM RR S 6IN SLVR	Silver	Grey	Grey	6	REAR

### DRIVE WHEEL INFORMATION

Drive Wheel Part #	Description	Hub Color	Tire color	Matching Frame Color	CHAIR LOCATION
125455	WHEEL DD 13 BLK 5 SPK GRY TIRE	Black	Grey	Black Frame	Center
126539	WHEEL DD 13 BLK 5 SPK BLK TIRE	Black	Black	Black Frame	Center
106468	DRIVE WHEEL ASSY 13	Silver	Grey	Silver Frame	Center