


QTS Tilt Mechanism Part History:


Generation I

Before 12/15/95



Generation II


12/15/95 to 11/15/96



Generation III


11/15/96 to present

(Dual Mech Lock option)



1 X 920411 *
Part includes:
Silver color
(2-Mech Locks)-342039
& Black color
(1- Gas Damper)-920404

(Gas Cylinder option)




2 X 920371
Black color
(Gas Strut)-342064

or

Generation IV (growing frame)


03/25/99 to present

(Dual Mech Lock option)



1 X 920950 *
Part includes:
Silver color
(2-Mech Locks)-342039
& Black color
(1- Gas Damper)-920404

(Gas Cylinder option)



2 X 920371
Black color
(Gas Strut)-342064

or

* This part # includes both Mech Locks and a Gas Damper when ordered, even though only one or the other is needed.