

SEAT FRAME - TILT & RECLINE

[10/2011]

NOTE: EFFECTIVE 4/21/03 STARTIING S/N S64-034729.

Pos.	Item Number	Description	Pos.	Item Number	Description
1	223121	3/8-16 x .75 HEX HD C/S W/PATCH	33	542506	20W TILT CHNNL WLDMT (.5) TUBE
2	262141	3/8 LOCK WASHER	33	542507	22W TILT CHNNL WLDMT (.5) TUBE
3	349598	BEARING IGUS 5/8 x 3/4 x 1/2	33	542508	24W TILT CHNNL WLDMT (.5) TUBE
4	481466	REC PIVOT BOSS	33A	542530	TILT CHANNEL WLDMNT 14W LONG
5	117829	REAR PIVOT POST ADJ	33A	542531	TILT CHANNEL WLDMNT 15W LONG
6	226064	1/4-28 x 3/4 BH C/S BLK W/PTCH	33A	542532	TILT CHANNEL WLDMNT 16W LONG
7	106943	WASHER FLAT AN 1/4 ROHS	33A	542533	TILT CHANNEL WLDMNT 17W LONG
8	227110	1/4-28 x 1 BH C/S BLK ZN	33A	542534	TILT CHANNEL WLDMNT 18W LONG
9	542162	REAR PIVOT POST ADJ MOUNT	33A	542535	TILT CHANNEL WLDMNT 19W LONG
10	481419	1/4-28 CHANNEL NUT 1.75" SP	33A	542536	TILT CHANNEL WLDMNT 20W LONG
11	221097	1/4-28 x 5/8 HEX C/S BZ	33A	542537	TILT CHANNEL WLDMNT 22W LONG
13	505929	HOME SWITCH BRACKET	33A	542538	TILT CHANNEL WLDMNT 24W LONG
14	481418	1/4-28 CHANNEL NUT .80" SP	34	567003	END PLUG OVAL
15	105787	PRSD A/A LATCH BLOCK NUT ASSM	35	139449	HOME SWITCH TILT ONLY
16	226218	5/16-18 x 1 3/4" BH C/S BLK	36	221090	1/4-28 x 3/8 BH C/S BLK
17	106944	WASHER FLAT THIN AN 1/4 ROHS	12,34	660180S	FLRD SEAT RAIL ASSM STD L S626
18	221012	1/4-28 x .44 SHCS BLK W/PL	12,34	660181S	FLRD SEAT RAIL ASSM STD R S626
19	226064	1/4-28 x 3/4 BH C/S BLK W/PTCH	12A,34	660182S	FLRD SEAT RAIL ASSM LNG L S626
20	202300	1/4-28 NYLOCK NUT THIN BLK ZN	12A,34	660183S	FLRD SEAT RAIL ASSM LNG R S626
21	106943	WASHER FLAT AN 1/4 ROHS	1-6,8-11,21	885115	REAR PIVOT POST ADJ ASSM TR
22	226247	1/4-28 x 5/8 B/H C/S BLACK G5	17-33	885220	TILT CHNNL WLDMNT ASSM 14W STD
23	221092	1/4-28 x 7/16 BH CS BLK W/PTCH	17-33	885221	TILT CHNNL WLDMNT ASSM 15W STD
24	740142	SUPPORT STRAP T & R ACTUATOR	17-33	885222	TILT CHNNL WLDMNT ASSM 16W STD
25	481106-QU	1/4-28 CHANNEL NUT	17-33	885223	TILT CHNNL WLDMNT ASSM 17W STD
26	222230	5/16-18 x 3/4 HEX C/S G8	17-33	885224	TILT CHNNL WLDMNT ASSM 18W STD
27	263100	5/16 AN WASHER FLAT	17-33	885225	TILT CHNNL WLDMNT ASSM 19W STD
28	202116	5/16-18 NYLOCK NUT THIN BLK	17-33	885226	TILT CHNNL WLDMNT ASSM 20W STD
29	202117	5/16-18 NYLOCK NUT THICK BLK	17-31,32A,33	885227	TILT CHNNL WLDMNT ASSM 22W STD
30	226188	5/16-18 x 1 FH C/S	17-31,32A,33	885228	TILT CHNNL WLDMNT ASSM 24W STD
31	502270	PIVOT BRKT ADJ REAR	17-32,33A	885229	TILT CHNNL WLDMNT ASSM 14W LNG
32	505609	REAR LOWER MULTI CLEVI WLDMNT	17-32,33A	885230	TILT CHNNL WLDMNT ASSM 15W LNG
32A	505623	REAR LWR MLT CLVS MNT TRAY LRG	17-32,33A	885231	TILT CHNNL WLDMNT ASSM 16W LNG
33	542500	14W TILT CHNNL WLDMT (.5) TUBE	17-32,33A	885232	TILT CHNNL WLDMNT ASSM 17W LNG
33	542501	15W TILT CHNNL WLDMT (.5) TUBE	17-32,33A	885233	TILT CHNNL WLDMNT ASSM 18W LNG
33	542502	16W TILT CHNNL WLDMT (.5) TUBE	17-32,33A	885234	TILT CHNNL WLDMNT ASSM 19W LNG
33	542503	17W TILT CHNNL WLDMT (.5) TUBE	17-32,33A	885235	TILT CHNNL WLDMNT ASSM 20W LNG
33	542504	18W TILT CHNNL WLDMT (.5) TUBE	17-31,32A,33A	885236	TILT CHNNL WLDMNT ASSM 22W LNG
33	542505	19W TILT CHNNL WLDMT (.5) TUBE	17-31,32A,33A	885237	TILT CHNNL WLDMNT ASSM 24W LNG