

Product Bulletin

 Power

 Manual Folding

 Manual Rigid/Sports

 Zippie Series

 Seating & Positioning

Torx Fasteners Update

Brent Hatch (May 24, 2006)

Detail:

Based on market feedback Sunrise Medical has removed Torx fasteners (bolts) from the following Quickie Adult Manual and Sports products:

- Quickie GT
- Quickie AllCourt (Aluminum & Titanium)
- Quickie MatchPoint (Aluminum & Titanium)

Torx fasteners have been replaced with metric Allen Head fasteners (see image below).

Solution:

- 1) All production wheelchairs listed above will use metric Allen Head fasteners starting on the implementation date listed below.
- 2) Corresponding spare parts pages have been updated. Note, replacement metric Allen Head fasteners can be used on product sold before the implementation date.
- 3) Tool Kit for the Quickie GT has been updated with new wrenches (Sunrise part # 050151)

Implementation Date: May 24, 2006

