

GUARDIAN TREK 3/TREK 4 PARTS SUPPLEMENT




Guardian Trek

IN THIS SECTION

Tiller/Midbase (Trek 3).....	1.1	Front Wheel Assembly (Trek 3).....	1.6
Tiller/Midbase (Trek 4).....	1.2	Front Wheel Assembly (Trek 4).....	1.7
Rear Assembly (Trek Drive Pak).....	1.3	Rear Wheel Assembly (Trek 3&4).....	1.7
Tote Wheels (Trek 3&4).....	1.3	Tiller Locking Assembly/Tote Handle (Trek 3).....	1.8
Seat (Trek 3&4).....	1.4	Tote Handle/Bumper (Trek 4).....	1.8
Seat post w/Pin (Trek 3&4).....	1.4	Basket w/Mounting Hardware (Trek 3&4).....	1.9
Armrests w/Knobs (Trek 3&4).....	1.5	Tiller Locking Knob Assembly (Trek 3&4).....	1.9
Battery Pack w/Charger (Trek 3&4).....	1.5	Keys (Trek 3&4).....	1.9
Front Steering/Tierod Assembly (Trek 4).....	1.6		

**TILLER/MIDBASE
(GUARDIAN TREK 3)**



TREK0X07

ITEM	NUMBER	DESCRIPTION	QTY
1	100417	TILLER MIDBASE ASSY RED	1

ITEM	NUMBER	DESCRIPTION	QTY
1A	100418	TILLER MIDBASE ASSY BLUE	1

* NOTE: TILLER/MIDBASE ASSEMBLY INCLUDES FRONT WHEEL/FORK, LOCKING KNOB, CLAMP ASSEMBLY FOR BASKET, AND KEYS.

**TILLER/MIDBASE
(GUARDIAN TREK 4)**
(Effective October 1 2004)




TREK4X03

ITEM	NUMBER	DESCRIPTION	QTY
1	100419	G4 TILLER AND BASE ASSM RED	1

ITEM	NUMBER	DESCRIPTION	QTY
1A	100420	G4 TILLER AND BASE ASSM BLUE	1

* NOTE: TILLER/MIDBASE ASSEMBLY INCLUDES FRONT WHEEL/FORK, LOCKING KNOB, CLAMP ASSEMBLY FOR BASKET, AND KEYS.

**REAR ASSEMBLY (TREK DRIVE PAK)
(GUARDIAN TREK 3/TREK 4)**



TREK0X04

ITEM	NUMBER	DESCRIPTION	QTY
1	492111	REAR ASSY DRIVE PAK RED	1

ITEM	NUMBER	DESCRIPTION	QTY
1A	492112	REAR ASSY DRIVE PAK BLUE	1

TOTE WHEELS (GUARDIAN TREK 3/TREK 4)



Trek0X19

ITEM	NUMBER	DESCRIPTION	QTY
1-4	880987	ANTI-TIP WHEEL ASSEMBLY (PR)	

**SEAT
(GUARDIAN TREK 3/TREK 4)**



TREK0X08

ITEM	NUMBER	DESCRIPTION	QTY
1	888928	GRDN TREK FOLDING SEAT	1

SEATPOST AND PIN (GUARDIAN TREK 3/TREK 4)



TREK0X03


ITEM	NUMBER	DESCRIPTION	QTY
1	500347	GRDN TREK SEAT POST	1

ITEM	NUMBER	DESCRIPTION	QTY
2	330071	ADJUSTMENT PIN M8 x 52.5	1

FLIP-DOWN ARMREST W/KNOBS

(GUARDIAN TREK 3/TREK 4)

(Effective July 1 2004)




TREK0X05

ITEM	NUMBER	DESCRIPTION	QTY
1,2	888902	ARM SET W/KNOBS	

ITEM	NUMBER	DESCRIPTION	QTY
2	349049	KNOB M8 x 1.25 x 18	2

BATTERY PACK W/ CHARGER

(GUARDIAN TREK 3/TREK 4)




TREK0X12

ITEM	NUMBER	DESCRIPTION	QTY
2	503026	BAT PCK ACCESS DOOR	1
1,2	136074	BATTERY PACK	1

ITEM	NUMBER	DESCRIPTION	QTY
3,4	888901	OFF BOARD CHARGER	1
4	134121	POWER CORD BATT CHARGER YAMAHA	1
5*	BAT12	BATTERY 12 VOLT	

*NOTE: TO REPLACE INTERNAL BATTERIES ORDER QTY 2 OF BAT 12.


FRONT STEERING/TIEROD ASSEMBLY (GUARDIAN TREK 4)



TREK4X01

ITEM	NUMBER	DESCRIPTION	QTY
1-7	888974	G4 FRONT STEERING ASSEMBLY	

FRONT WHEEL ASSEMBLY (GUARDIAN TREK 3)




TREK0X10

ITEM	NUMBER	DESCRIPTION	QTY
1	780012	GRDN TREK 3 FRONT TIRE/RIM	1

ITEM	NUMBER	DESCRIPTION	QTY
1-7	888926	GRDN TREK 3 FRONT WHEEL ASSM	

FRONT WHEEL ASSEMBLY (GUARDIAN TREK 4)




TREK4X02

ITEM	NUMBER	DESCRIPTION	QTY
1	780013	GRDN TREK 4 FRONT TIRE/RIM	1

ITEM	NUMBER	DESCRIPTION	QTY
1-5	888973	G4 FRONT WHEEL W/HARDWARE KIT	

REAR WHEEL ASSEMBLY (GUARDIAN TREK 3/TREK 4)




TREK0X09

ITEM	NUMBER	DESCRIPTION	QTY
1	760001	GRDN TREK REAR TIRE/RIM	1

ITEM	NUMBER	DESCRIPTION	QTY
1-6	888927	GRDN TREK 3 REAR WHEEL ASSM	1

* NOTE: WHEN ORDERING REAR TIRE ASSEMBLIES, AS A PAIR, THEY MUST BE ORDERED AS A QUANTITY OF 2 (888927).


TILLER LOCKING ASSEMBLY/TOTE HANDLE (GUARDIAN TREK 3)



TREK0X16

ITEM	NUMBER	DESCRIPTION	QTY
1-4	888872	TILLER LOCKING ASSM	1


TOTE HANDLE/BUMPER (GUARDIAN TREK 4)



TREK4X04

ITEM	NUMBER	DESCRIPTION	QTY
1-5	888975	G4 TOTE HNDL/FRONT BUMPER ASSM	1


BASKET W/MOUNTING HARDWARE (GUARDIAN TREK 3/TREK 4)



TREK0X11

ITEM	NUMBER	DESCRIPTION	QTY
4	888923	GRDN TREK 3 BASKET W/O HDWR	


TILLER LOCKING KNOB ASSEMBLY (GUARDIAN TREK 3/TREK 4)



TREK0X15

ITEM	NUMBER	DESCRIPTION	QTY
1-3	888932	TILLER LOCKING KNOB ASSM	1

KEYS (GUARDIAN TREK 3/TREK 4)



TREK0X13

ITEM	NUMBER	DESCRIPTION	QTY
1	331090	KEYS	1
1	100424	KEYS (AFTER 8/1/05)	1

Note: Order P/N 331090 for scooters sold before 8/1/05 and use P/N 100424 for scooters sold after 8/1/05.