

LEFT FOOT RELEASE TILT RELEASE				
CHAIR WIDTH	DEPTH RANGE (ROLLER TO FOOT PEDAL DISTANCE)	LEFT CABLE LENGTH		
		TILT REL FOOT L SHRT 160032	TILT REL FOOT L MED 160033	TILT REL FOOT L LNG 160899
10	12-14	X	-	-
	15-16	X	-	-
	17-18	-	X	-
11	12-14	X	-	-
	15-16	X	-	-
	17-18	-	X	-
12	12-14	X	-	-
	15-16	-	X	-
	17-18	-	X	-
13	12-14	-	X	-
	15-16	-	X	-
	17-18	-	-	X
14	12-14	-	X	-
	15-16	-	X	-
	17-18	-	-	X
15	12-14	-	X	-
	15-16	-	X	-
	17-18	-	-	X
16	12-14	-	X	-
	15-16	-	X	-
	17-18	-	-	X
17	12-14	-	X	-
	15-16	-	-	X
	17-18	-	-	X
18	12-14	-	-	X
	15-16	-	-	X
	17-18	-	-	X

RIGHT FOOT RELEASE TILT RELEASE			
CHAIR WIDTH	DEPTH RANGE (ROLLER TO FOOT PEDAL DISTANCE)	RIGHT CABLE LENGTH	
		TILT REL FOOT R SHRT 160031	TILT REL FOOT R LNG 160898
10	12-14	X	-
	15-16	-	X
	17-18	-	X
11	12-14	X	-
	15-16	-	X
	17-18	-	X
12	12-14	X	-
	15-16	-	X
	17-18	-	X
13	12-14	X	-
	15-16	-	X
	17-18	-	X
14	12-14	X	-
	15-16	-	X
	17-18	-	X
15	12-14	X	-
	15-16	-	X
	17-18	-	X
16	12-14	X	-
	15-16	-	X
	17-18	-	X
17	12-14	X	-
	15-16	-	X
	17-18	-	X
18	12-14	X	-
	15-16	-	X
	17-18	-	X

COMBINATIONS:	DESCRIPTION:	LEFT CABLE	RIGHT CABLE
	LFT SHRT RHT SHRT	160032	160031
	LFT SHRT RHT LNG	160032	160898
	LFT MED RHT LNG	160033	160898
	LFT MED RHT SHRT	160033	160031
	LFT LNG RHT LNG	160899	160898
	LFT LNG RHT SHRT	160899	160031

REAR WHEEL SIZE	ROLLER POS	FOOT PEDAL POS	DEPTH RANGE (ROLLER TO FOOT PEDAL DISTANCE)
12-16	1	14	13
		15	14
		16	15
		17	16
		17	16
	2	15	13
		16	14
		17	15
		19	17
		19	17
	3	15	12
		16	13
		17	14
		19	16
		19	16
4	16	12	
	17	13	
	19	15	
	19	15	
	19	15	
20-24	1	17	16
		18	17
		19	18
	2	18	16
		19	17
		19	17
	3	18	15
		19	16
	4	19	15
		19	15

