

PULL WHEEL LOCK AVAILABILITY (Offset, Curved, and Straight Axle Plates)						
Frame		Small	Medium	Long	X Long	
Rear Wheel Size	Rear Seat To Floor Height					
20	13.0	Forward Mount Only				
	13.5					
	14.0					
	14.5					
	15.0					
	15.5					
	16.0					
	16.5					
	17.0			Reverse Mount Only		
	17.5					
18.0						
18.5						
22	13.5	Forward Mount Only				
	14.0					
	14.5					
	15.0					
	15.5					
	16.0					
	16.5					
	17.0					
	17.5					
	18.0					
18.5						
19.0						
19.5						
20.0						
24/25	14.5	N/A	Forward Mount Only			
	15.0					
	15.5					
	16.0					
	16.5					
	17.0					
	17.5					
	18.0					
	18.5					
	19.0					
19.5						
20.0						
20.5	Forward Mount Only					
21.0						
21.5						
26	15.5	N/A	N/A	Forward Mount Only		
	16.0					
	16.5					
	17.0					
	17.5					
	18.0					
	18.5					
	19.0					
	19.5					
	20.0					
20.5	Forward Mount Only					
21.0						
21.5						
22.0						

PULL WHEEL LOCK AVAILABILITY (Amputee Axle Plate)							
Frame		Small	Medium	Long	X Long		
Rear Wheel Size	Rear Seat To Floor Height						
20	14.0	Forward Mount Only **** N/A with SPHA & SPHA Flip Back Arms	Forward Mount Only **** N/A with SPHA & SPHA Flip Back Arms	Reverse Mount Only	Reverse Mount Only		
	14.5						
	15.0						
	15.5						
	16.0						
	16.5						
	17.0					Forward Mount Only	Forward Mount Only
	17.5					N/A with SPHA & SPHA Flip Back Arms	N/A with SPHA & SPHA Flip Back Arms
	18.0						
	18.5						
19.0							
22	14.5	N/A	N/A	Reverse Mount Only	Reverse Mount Only		
	15.0						
	15.5						
	16.0						
	16.5						
	17.0						
	17.5						
	18.0						
	18.5						
	19.0						
19.5	Forward Mount Only	Forward Mount Only					
20.0	N/A with SPHA & SPHA Flip Back Arms	N/A with SPHA & SPHA Flip Back Arms					
24/25	15.5	N/A	N/A				
	16.0						
	16.5						
	17.0						
	17.5						
	18.0						
	18.5						
	19.0						
	19.5						
	20.0						
20.5			Reverse Mount Only	Reverse Mount Only			
21.0							
26	17.0	N/A	N/A				
	17.5						
	18.0						
	18.5						
	19.0						
	19.5						
	20.0						
	20.5						
	21.0						
	21.5						
22.0							